

## LAB SCALE EXTRUSION ALTERNATIVE PROTEINS AND MEAT ANALOGUES

Brabender becomes  
part of Anton Paar.  
Two pioneers,  
one mission.

## FACTS & FIGURES



ESTABLISHED IN  
**1922**



HEADQUARTERS  
**IN GRAZ / AUSTRIA**



**4,400+**  
EMPLOYEES



OWNED BY THE CHARITABLE  
**SANTNER FOUNDATION**



**14.5 % INVESTMENT**  
IN RESEARCH AND DEVELOPMENT  
FROM ANNUAL TURNOVER  
ANTON PAAR GMBH




**ALL CRITICAL COMPONENTS**  
MANUFACTURED IN-HOUSE

BRABENDER HAS BEEN INTEGRATED  
AS A BRAND OF ANTON PAAR



**Anton Paar**

**Brabender®**  
A brand of **Anton Paar**



## **BRABENDER PRODUCT LINE**

---

- › SERVING INDUSTRY & RESEARCH SINCE 1923
- › PART OF THE ANTON PAAR GROUP SINCE AUGUST 2023

## **PRODUCT PORTFOLIO**

- › MODULAR TORQUE RHEOMETERS
- › MEASURING MIXERS
- › FLOUR, DOUGH & STARCH ANALYZERS
- › VISCOMETERS
- › MOISTURE, DENSITY, OIL ABSORPTION & RELAXATION ANALYZERS
- › LAB-SCALE EXTRUDERS & PERIPHERY
- › LAB MILLS
- › PROCESS SENSORS
- › LABORATORY SOFTWARE

# ... AND IS NOW SERVING INDUSTRY AND RESEARCH WITH NUMEROUS APPLICATIONS



**FOOD & FEED**



**PLASTICS & RUBBER**



**OTHERS**

# THE BRABENDER HISTORY




**1923**

Company foundation by Carl Wilhelm Brabender. Corporate purpose: **Major repairs of electric motors and transformers.**




**1928**

Invention of the **Farinograph** as world's first torque measuring instrument for determining the **quality of flour.**



**1936**

Development of the **Plastograph**, based of the Farinograph's measuring principle. First measuring mixer on the market for determining **plastics and rubber quality.**



**1965**

Brabender enters the extrusion market with **single screw extruders for realistic process simulation at laboratory scale.**



**1970**

Carl Wilhelm Brabender's wife **Martha Brabender** takes over the company after the death of her husband.

# THE BRABENDER HISTORY



**2001**

Brabender develops the market's **first twin-screw extruder with hinged and openable liner**, optimized for process monitoring and cleaning.



**2012**

Brabender's first rheological instrument for **rapid quality testing: The GlutoPeak** was developed to measure gluten quality.



**2015**

First devices equipped with **MetaBridge software**, allowing customers to retrieve measurement readings from any device and location.



**2018**

ViscoQuick – **first universal viscometer** developed by Brabender. Measures starch-based products and fluids with different viscosities.




**2023**

Brabender celebrates its **100th anniversary** and is integrated into the **Anton Paar Group**





OUR PRODUCTION FACILITY IN DUISBURG, GERMANY



## GLOBAL BRABENDER LOCATIONS

# BRABENDER PRODUCT LINE FOOD & FEED



## Flour & Dough Analyzers

ICC, AACCI, ISO standard methods:

- › **FarinoGraph**: Flour water absorption and kneading characteristics of dough
- › **ExtensoGraph**: Dough elasticity and processing properties
- › **Amylograph-E**: Enzyme activity and baking properties of flour

Rapid method:

- › **GlutoPeak**: Aggregation properties and strength of gluten in flour (rapid method)



## Starch & Gluten Analyzers

Starch analyzers:

- › **Viscograph-E**: Starch gelatinization and enzyme activity (ICC, AACCI, standard method)
- › **ViscoQuick**: Starch gelatinization and enzyme activity (rapid method)

Gluten analyzers:

- › **GlutoPeak**: Vital gluten aggregation properties and strength (rapid method)
- › **Glutograph-E**: Gluten elasticity



## Mills & Moisture Analyzers

Equipment for sample preparation

**Lab Mills:**

Break mills and separators at laboratory scale for preparing samples for further analyses

**Moisture analyzer:**

- › **MT-CA**: Automated oven drying method with integrated balance and capacity for testing 10 samples at once (based on ICC and ISO standard method)

# BRABENDER PRODUCT LINE FOOD



## Lab-scale Extruders

- › Stand-alone single-screw extruders
- › Stand-alone twin-screw extruders
- › Single and twin-screw extruder attachments for modular torque rheometers
- › Application-specific dies and screws
- › Downstream equipment for conveying, cooling and cutting of extrudates

Applicable for e. g. snackfood, pasta, flatbread, meat analogues, feed pellets



## Laboratory Software

**MetaBridge:** Browser-based standard instrument software

- › Touch-screen optimized operation
- › Compatible with all desktop or mobile terminal devices
- › Compliance with standard and individually defined methods
- › Comparisons via reference and correlation add-on functions
- › Measured values transfer to other Brabender devices and 3<sup>rd</sup> party systems

# SINGLE SCREW AND TWIN SCREW EXTRUDERS AND ACCESSORIES

## Single screw extruders

Complete systems and attachments for the rheological drive units with 19 and 30 mm screw diameter



## Twin-screw extruders TwinLab Series

12/36 attachment  
20/40 system and attachment  
30/40 system



# SINGLE VS. TWIN SCREW EXTRUDERS

## **When to use a single-screw extruder?**

- › Customer only has one material or premix/compound
- › Customer does not want to actively adapt the recipe during the process
- › Customer requires a higher maximum process pressure (700 bar)
- › Save money with simplified setups


## **When to use a twin-screw extruder?**

- › Required if you want to mix more than 2 material
- › Throughput greater than 20 kg/h required
- › Future-proof and flexible solution when money is no issue



# SINGLE VS. TWIN SCREW EXTRUDERS

Degree of fill for single- and twin-screw




Schematic!

# CLAM SHELL DESIGN

## User-friendly design of Brabender twin-screw extruders

- › Openable liner
- › Easy cleaning
- › Process transparency





## Liner opening for screw design optimization

- **Opening of the liner enables the operator to have a look insight the process.**
- **Further benefit: cleaning process much easier**
- Identification of fully and half-filled zones
- Influence of temperature on the raw material
- Determination of the degree of mixing


⇒ **Optimization of screw design**



## Impact of screw elements on raw material




Mixture/ cold forming



Gelatinization/ plastification




Starch destruction



Co-rotating,  
intermeshing,  
twin-screw

# Extrusion processing... what about meat analogues?



CCO public domain/ pixabay




# Most meat analogues are extruded products...




# Extrusion process - HME

## meat analogue – extrusion process & modular cooling die






## lamination



# Extrusion process - HME


## meat analogues – HME

Starting procedure	Parameter determination	Sample production	Stopping procedure
Low shear + heat	Increase: shear + heat stepwise	Optimal conditions	Low shear + heat
	Under cooked		
	Over cooked		




$$\text{Relative fibre length} = \frac{L}{W}$$

$L = \text{Fibre length}$

$W = \text{Fibre width}$



## meat analogues – LME

Starting procedure	Parameter determination	Sample production	Stopping procedure
Low shear + heat	Increase: shear + heat stepwise	Optimal conditions	Low shear + heat
	Under cooked		
	Over cooked		



Soy chunks

# Set up & equipment

Marel demo center facilities:





# Set up & equipment

**Brabender extrusion line:**




**Twin screw extruder TwinLab-F 20/40  
+ application-specific screws**



**Round dies for LME**



**Cutting device for LME**



**Modular cooling die for HME**



**Conveyor belt for HME**



# Formulations & products

- **LME:**
- soy protein concentrate: „alpha-8“ (Solae), 67% protein content
  - water
- **HME:**
- soy protein concentrate: „alpha-8“ (Solae), 67% protein content
  - water
  - **Burger patty mix: „vegan burger mix“ (DP&S):**
    - emulsion
    - binding system
    - spice mixture
  - **Schnitzel mix: „vegan chicken schnitzel“ (DP&S):**
    - emulsions
    - binding system
    - spice mixture
    - coating system

# Processes:

## burger patties, chicken-type burgers, schnitzels & chicken-type pieces



extrusion line



LME extruded  
products  
of different size



HME  
extrusion

# Process: Application 1 burger patties – tasting

Comparison of the marel-Brabender-products to commercially available products:



→marel-Brabender patty was smooth texture in comparison to other products  
→marel-Brabender patty and Garden Gourmet patty showed best results at the tasting

# Process: Application 2 chicken-type burgers, schnitzels & chicken-type pieces



fried schnitzels



fried chicken-type filet products

# Research & Applications

## Pulse flour Extrusion



RQ: Is it possible to produce a whole BURGER made from HME stripes?


### Setup:

- twinscrew extruder
- modular cooling die



Raw materials:	Protein content	Carbo-hydrates	Fibre	Fat
Soy protein concentrate (supplier: Solae)	67 %	data not available	data not available	8 %
SMART PRO Faba bean (supplier: Müller's Mühle)	58,3 %	7,4 %	14,3 %	4,6 %
SMART PRO Yellow pea (supplier: Müller's Mühle)	48,4 %	17 %	16,2 %	4,3 %

water absorption of different pulse flours analyzed in advance



### Extrusion process (HME; faba bean):

- total moisture content: 50%
- temperature profile: 40-70-80-110 °C
- cooling die temperature: 20°C
- screw speed: 100 rpm
- throughput: 2Kg/h


### Results:

- pure faba bean → low texturization, 'beany' taste
- mixtures with good texturization and sensory properties:
  - faba bean + y. pea
  - pea + oat





# Selection of pulses for HME and TVP extrusion



## Chickpea:

- No texturization and/or expansion
- due to low protein and high-fat content

## Red Lentils:

- No texturization through HME
- TVP performed well

## Yellow Pea:

- Suitable for HME

## Fava Bean:

- Good behavior for HME
- The best results due to the formation of texturized proteins
- Pure: beany taste

# Extrusion of pulse-based pasta




## SMART PRO Chickpea

Protein content	Carbohydrates	Fibre	Fat
36,8 %	24,5 %	15,0 %	11,7 %

### Setup:

- twinscrew extruder
- round die
- extrusion process:
  - temperature profile [°C]:  
HZ1 – HZ5: 50-80-100-120-130
  - total moisture content [%]: 25
  - throughput [kg/h]: 2
  - screw speed [U/min]: 700



## Results:

extruded



dried



cooked



Further optimization: using a degassing unit

# Sidestream utilization

## Stale bread



The amount of bread thrown away is far too high.

## Brewer's spent grains



Brewer's spent grain is a byproduct of the brewing industry. Due to its high protein content, it is often used for animal feed.

## Olive Stone Powder



Many olive are needed for nutrition and producing of olive oil. But the olive stones are discarded



# Sidestream utilization

## Setup

- Specially configured screw for High Moisture Extrusion
  - High shear
  - Different kneading blocks
  - Tooth element
  - Reconveying



# Sidestream utilization

## Setup

Parameter	SPC [%]	Side stream [%]	Moisture content [%]
Breadcrumbs	70	30	8,04
Olive stone powder	70	30	7,85
Brewer's spent grains	70	30	27,2
Soy protein concentrate	100	0	7,44



# Sidestream utilization

Soy protein  
concentrate

Stale bread


Brewer's  
spent grains

Olive Stone Powder



# AquaFeed extrusion

- Optimization of fish feed formulations
- Substitution of fish meal by insect meal
- Increasing fat content by
  - changing process parameters
  - adding different silica types
- Exhibition + poster presented at...
  - Aquaculture Europe 2021 (Portugal)
  - VICTAM 2022 (Utrecht)
- Presentation at AquaFeed Extrusion Conference 2022 (Utrecht)
- Silica project to be continued
- Exhibition + poster (+ presentation) at Aquaculture Europe 2022 (Italy)



# Confectionary


- It's hard to imagine confectionery shelves without the very common fruit laces,
- which come in a wide variety of shapes and colors with various surface treatments to be offered.
- From waxed, or sprinkled with sugar, to a sour sprinkling with a sugar-acid mixture, a lot is possible.
- Products filled with sugar pastes are just as varied on offer as products made from liquorice, which often also come with a sugar coated surface
- No matter how great the variety of products, there is one manufacturing process behind it: extrusion.





# Confectionary

- Challenges:
  - Dosing of high viscous slurry
  - Dosing low amounts of colour and flavor




Rotary Piston Pump for high viskos  
mass (up to 2 million mPas)



# Confectionary

## Co-Rotating vs. Counter-Rotating

- Benefit Co-Rotating:
  - very flexible regarding the application
  - 2 uniform screws
  - High shear possible
- Benefit Counter-Rotating:
  - Very low shear possible



Arlington Machinery ([www.arlingtonmachinery.com](http://www.arlingtonmachinery.com)), 29.06.2023

# Confectionary

- Complete Setup with
  - Slurry Pump (Rotary Piston Pump)
  - Colordosing
  - Flavordosing
  - Degassing Station
  - Round Die Head





# Confectionary

## Process Parameter

- Temperature Profil:
  - HZ1 – HZ5: 90-135-110-90-90
- Total throughput:
  - 4 – 5 kg/ h

Screwspeed:

- 100 rpm

Vakuum:

- -0,35 bar

Schmelztemperatur:

- Approx. 110 °C



# REFERENCES R&D INSTITUTIONS

<p>Hochschule Düsseldorf University of Applied Sciences</p> 	<p>Hochschule Fulda University of Applied Sciences</p> 			
 <p>versuchsanstalt für getreideverarbeitung</p>	 <p>UNIVERSITÀ DEGLI STUDI DI MILANO</p>	 <p><b>AIMPLAS</b> PLASTICS TECHNOLOGY CENTRE</p>		
<p>WASHINGTON STATE UNIVERSITY</p> 	<p><b>NDSU</b> NORTH DAKOTA STATE UNIVERSITY</p>	 <p><b>RISE</b>   Research Institutes of Sweden</p>	 <p><b>WAGENINGEN</b> UNIVERSITY &amp; RESEARCH</p>	
	 <p><b>ucc</b> Coláiste na hOllscoile Corcaigh, Éire University College Cork, Ireland</p>		 <p>UNIVERSIDAD NACIONAL DE COLOMBIA</p>	

# REFERENCES FOOD



# REFERENCES FOOD



# APPLICATION LABS

**Get support from experts in quality control and R&D**



**Run trials with equipment for rheological and other measurements, extrusion and more**





# APPLICATION LABS – remote lab

- Lab tours
- Device demos
- Material analysis
- Extrusion trials
- Product & process development



## Benefits

- No travel restrictions
- Limited costs
- Customers from different subsidiaries can attend
- Recording
- Live interaction with applic. team



# CUSTOMIZATION

## Custom-tailored solutions

- › Experts in technology and application
- › Customized components and systems





**Get in touch**

[www.anton-paar.com](http://www.anton-paar.com)

# EXTRUDER DIE HEADS AND SCREWS

- › Application-specific dies for the production of extrudates with different shapes
- › Screws for single screw extruder tailored to specific applications
- › Modular screws for twin screw extruders



# EXTRUDER DIE HEADS –EXAMPLES

- › Round die head
- › Ribbon die head
- › Multi-strand die head
- › Modular cooling die

